

JATI AUTOMASI ENTERPRISE

(Co. no 001773550-P)

No. 9B Jalan SG ¾, Taman Sri Gombak, 68100 Batu Caves, Selangor Darul Ehsan.
Tel : 03-61859326 Fax : 03-61871294



DIATERMO 106 monopolar – 50 watt

Technical Specifications:

Output power:
CUT 50W/400 ohm
COAG 40W/400 ohm
BLEND 45W/400 ohm
Working frequency: 600 KHz +/- 10%
Size: 180 x 180 x h 100 mm
Weight: 2.8 kg
Operating Voltage: 115/230 V - 50/60 Hz switchable
Electric input power: 90 VA
Fuse: 2 x 1 AT
Alarm for overheating
Footswitch operating
Norms: EN60601-1 Class I CF, EN 60601-1-2 93/42/EEC (Class IIb)
Warranty: 12 months



Description:

These HF electrosurgery units, are particularly suitable for small operations and are recommended for veterinary, dentistry, first aid and dermatology work. They offer precision of cut and immediate coagulation and give excellent results in a large range of applications: from uterine coagulation to microcoagulation in dermatology. Moreover you can choose from a wide variety of electrodes for any need. The modes of operation, CUT, COAG BLEND and MICRO (30517 only), can be easily selected by a switch on the front panel. The complete isolation between patient and the unit and strict respect of international norms guarantee total safety.

Standard accessories:

KIRKY single-use handle without buttons
Operative foot-switch
Flat neutral plate with cable (30563)
Kit of 6 short autoclavable electrodes, Ø 2.4 mm, length 5 cm
Power cable 2 m - Shuko plug
User manual: available in English, Italian, Spanish, French