

# JATI AUTOMASI ENTERPRISE

(Co. no 001773550-P)

No. 9B Jalan SG ¾ , Taman Sri Gombak, 68100 Batu Caves, Selangor Darul Ehsan.  
Tel : 03-61859326 Fax : 03-61871294

## HYDRAULIC DELIVERY BED – LIGHT VERSION MODEL DB 2020 (LV)

### Product Features:

- Convertible from bed to delivery table
- Steel construction finishing with scratch and chip resistant of hard epoxy powder coating
- Manual ratchet adjustable backrest
- Hydraulic system height adjustable foot pedal at both side of bed
- Instant trendelenburg & reverse trendelenburg position assisted by gas spring with a control lever on both side of bed
- Retractable leg section - slides smoothly on steel bearings under the sectional base and locks automatically
- Two hinged section PVC padded foam mattress of 75mm thick & a 100mm thick for leg section
- Inclusive of IV pole, a pair of side rails, a detachable head board, a stainless steel bowl & a pair of lithotomy stirrups and straps
- Mobile on 4 x 125mm swivel castors with 2 brakes diagonally fitted
- Removable foot step

### Technical Specification:

- Hi-Lo: 610 - 950(H)
- Overall size: 1200(L) stowed / 1930(L) extended 880(W)

