

# JATI AUTOMASI ENTERPRISE

(Co. no 001773550-P)

No. 9B Jalan SG ¾, Taman Sri Gombak, 68100 Batu Caves, Selangor Darul Ehsan.  
Tel : 03-61859326 Fax : 03-61871294

## **ELECTRIC DELIVERY BED MODEL DB 4000**

### **Product Features:**

- Convertible from bed to delivery unit
- With easy detachable PE blow head & leg boards
- Hi-Lo, backrest, leg rest and trendelenburg position controlled by electrical motorized
- Foot pedal controls unit to allow a hands free for medical personnel
- Retractable leg section
- Manual Override CPR release lever at both sides for backrest
- Mobile on 4 x 150mm dia. Tente (Germany) central locking system
- Removable molded mattress base
- Double PE blow side rail with gas spring assisted adjustment
- Easy fold-away side rails
- Delivery position, Squatting position, spreader footrest position, lithotomic position

### **Electrical System:**

- LINAK DC-actuator system (TUV & UL approved & EC60601-1)
- 24VDC out-put
- Built in back-up battery

### **Technical Specification:**

- Height adjustment: 560 - 940mm
- Backrest : 0o - 80o
- Trendelenburg: 0o - 8o
- Overall size(mm): 1448(L) stowed/2108(L) extended x 813(W) x 600-990(H)
- Safety load capacity: 230kg UDL

# JATI AUTOMASI ENTERPRISE

(Co. no 001773550-P)

No. 9B Jalan SG ¾, Taman Sri Gombak, 68100 Batu Caves, Selangor Darul Ehsan.

Tel : 03-61859326 Fax : 03-61871294



## KIT STAND

- \* IV drip stand
- \* 2 Swivel padded leg support
- \* 2 Spreaded footrest
- \* 1 Patient step
- \* 2 Squatting support
- \* 1 Stainless steel fluid collection bowl