

# JATI AUTOMASI ENTERPRISE

(Co. no 001773550-P)

No. 555 C Jalan Samudra Utara 1, Taman Samudra, 68100 Batu Caves, Selangor Darul Ehsan.  
Tel : 03-61859326 Fax : 03-61871294

## EXAMINATION LIGHT - ANGLE POISE MODEL WM 700

### Features / Specification:

- With spring adjustable tension angle poise arm & lamp shade
- C/w a 40watt bulb.
- Mobile on 4 x 50mm castors

### Make:

- Constructed of tubular steel leg frame coated with hard epoxy powder coating

