

JATI AUTOMASI ENTERPRISE

(Co. no 001773550-P)

No. 555 C Jalan Samudra Utara 1, Taman Samudra, 68100 Batu Caves, Selangor Darul Ehsan.
Tel : 03-61859326 Fax : 03-61871294

EXAMINATION TABLE FOR GYNAE MODEL EC 5000

Features:

- All parts are made by tooling to ensure uniform & outstanding quality
- Fixed Height
- 75mm thick PVC mattress
- Ratchet backrest
- Come with a pair of lithotomy stirrups and straps

Make:

- Made of tubular steel frame work fitted on rubber studs
- Finishing with a layer of chip and scratch resistant baked epoxy coat. (cream colour)

Size

- 1830mmL x 610mmW x 760mmH

*Available in Complete with a stainless steel basin. (EC 5001)

