

JATI AUTOMASI ENTERPRISE

(Co. no 001773550-P)

No. 9B Jalan SG ¾, Taman Sri Gombak, 68100 Batu Caves, Selangor Darul Ehsan.
Tel : 03-61859326 Fax : 03-61871294

HOSPITAL BEDSIDE RACK MODEL WP 1001

Product Features

- Steel construction
- Scratch & chip resistant epoxy coating
- Chipboard laminated surface with all round T-Frame rubber trimmings edge table top
- 50mm swivel castors
- Dim (mm): 490(W) x 406(D) x 815(H)

