

JATI AUTOMASI ENTERPRISE

(Co. no 001773550-P)

No. 555 C Jalan Samudra Utara 1, Taman Samudra, 68100 Batu Caves, Selangor Darul Ehsan.
Tel : 03-61859326 Fax : 03-61871294

TABLE, INSTRUMENT S/S MAYO MODEL CT 4000 MODEL CT 4010

Features:

- Tubular steel hollow section pipe with chrome support rim for an s/steel tray
- With adjustable height by releasing & screw knob to lock.
- Mobile on 50mm castors.

Make:

- Made of tubular hollow section pipe.
- Finishing of hard epoxy powder coating cream color.

Tray size:

- 490mmL x 320mmW (CT 4000)
- 542mmL x 418mmW (CT 4010)

Overall Height:

- 810 – 1230mmH

