

JATI AUTOMASI ENTERPRISE

(Co. no 001773550-P)

No. 555 C Jalan Samudra Utara 1, Taman Samudra, 68100 Batu Caves, Selangor Darul Ehsan.
Tel : 03-61859326 Fax : 03-61871294

MEDICAL NURSERY CRIB MODEL DC1010

Product Features

- ✚ Steel Frame
- ✚ Scratch and chip resistant epoxy coating
- ✚ Clamp lock tilting device
- ✚ Closed cabinet
- ✚ 50mm thick PVC mattress
- ✚ Mobile on 75mm swivel castors
- ✚ Dim(mm): 765(L) x 445(W) x 865(H)



MEDICAL NURSERY CRIB MODEL DC1000

Product Features

- ✚ Steel Frame
- ✚ Scratch and chip resistant epoxy coating
- ✚ Position Ratchet
- ✚ Open shelf
- ✚ Mobile on 50mm thick PVC mattress
- ✚ 75mm swivel castors
- ✚ Dim (mm): 785(L) x 405(W) x 890(H)



MEDICAL NURSERY CRIB MODEL DC 1011

Product Features

- ✚ Steel Frame
- ✚ Scratch and chip resistant epoxy coating
- ✚ Position Ratchet
- ✚ Closed cabinet & open shelf
- ✚ 50mm thick PVC mattress
- ✚ Mobile on 75mm swivel castors
- ✚ Dim (mm): 785(L) x 405(W) x 890(H)

