

JATI AUTOMASI ENTERPRISE

(Co. no 001773550-P)

No. 9B Jalan SG ¾, Taman Sri Gombak, 68100 Batu Caves, Selangor Darul Ehsan.
Tel : 03-61859326 Fax : 03-61871294

MEDICAL NURSERY RADIANT HEATER MODEL DC 4000

Features: -

- Steel hollow section pipe frame.
- Radiant box of an Infrared radiant ceramic heater
- Height adjustment by gas spring tension & lock by a cap screw knobs.
- Individual On/Off switch with indicator.

Castor: -

- Four 50mm swivel castors

Finishing: -

- Hard epoxy powder coating – grey colour

Overall Size: -

- 400mmL x 700mmW x (1320 - 1700mm)H

Base Dim: -

- 400mmL x 700mmW

