

JATI AUTOMASI ENTERPRISE

(Co. no 001773550-P)

No. 555 C Jalan Samudra Utara 1, Taman Samudra, 68100 Batu Caves, Selangor Darul Ehsan.
Tel : 03-61859326 Fax : 03-61871294

MEDICAL IV TRANSFUSION STAND MODEL WM 200

- Tubular steel frame stand with height adjustable & screw knob to lock.
- C/w a double hook pole - chrome plated
- Base of PU 5 legged mobile on 50mm swivel castors.

