

JATI AUTOMASI ENTERPRISE

(Co. no 001773550-P)

No. 9B Jalan SG ¾, Taman Sri Gombak, 68100 Batu Caves, Selangor Darul Ehsan.
Tel : 03-61859326 Fax : 03-61871294

OVERBED TABLE - ADJUSTABLE HEIGHT MODEL WP 2100

Product Features

- Steel construction
- Scratch & chip resistant epoxy coating
- Knob locking height adjustable by classic tension spring (790mm to 1080mm)
- Chipboard laminated surface with all round T-Frame rubber trimmings edge table top
- 50mm swivel castors
- Table Top Dim(mm): 810(L) x 360(W)
- Dim(mm) : 825(L) x 406(W) x 760-1080(H)



OVERBED TABLE – ADJUSTABLE HEIGHT MODEL WP 2200

Product Features

- Steel construction
- Scratch & chip resistant epoxy coating
- Fingertip Height Adjustable by roller tension spring (760mm to 1080mm)
- Chipboard laminated surface with all round T-Frame rubber trimmings edge table top
- 50mm swivel castors
- Table Top Dim(mm): 810(L) x 360(W)
- Dim(mm) : 825(L) x 406(W) x 760-1080(H)

