

JATI AUTOMASI ENTERPRISE

(Co. no 001773550-P)

No. 555 C Jalan Samudra Utara 1, Taman Samudra, 68100 Batu Caves, Selangor Darul Ehsan.
Tel : 03-61859326 Fax : 03-61871294

OMRON NE-U17 MESH NEBULIZER

Features

- *High volume nebulizer for humidification or medication therapy
- *Variable output up to 4ml/min
- *For clinical applications, capable of running continuously for up to 72 hours.
- *Very quiet
- *Backlit display shows settings and timer

Specifications

- *Power Source: 230 VAC, 50 Hz
 - *Power Consumption: Approximately 80 VA
 - *Fuse: T2A L250V
 - *Ultrasonic Frequency: Approximately 1.7 MHz
 - *Nebulization Rate*: Approximately 0 to 3ml (adjustable)
 - *Particle Size*: 4.4 μ MMAD, based on cascade impaction particle size measurement method (according to EN13544-1 European nebuliser standard)
 - *Air Volume: Maximum of 17 l/min
 - *Amount of Cooling Water: Approximately 375 ml
 - *Capacity of Medication Cup: Approximately 150 ml (min.5 ml)
 - *External Dimensions: Approximately 276 (w) x 243 (H) x 226 (D)mm
 - *Weight of the Main Unit: Approximately 4.0 kg
 - *Operating Conditions: 5C to 40C (41F to 73.4F), 30 to 85% RH
 - *Accessories Included: Inhalation hose M (70cm with a cuff), mouthpiece set, 2 medication cups, power cord, instruction manual (with warranty card)
- *May change depending on the type of medication
*Specifications and appearance is subject to change without prior notice
*Warranty:1 Year

