

JATI AUTOMASI ENTERPRISE

(Co. no 001773550-P)

No. 9B Jalan SG ¼ , Taman Sri Gombak, 68100 Batu Caves, Selangor Darul Ehsan.
Tel : 03-61859326 Fax : 03-61871294

Omron NE-C801 Compressor Nebulizer

Features

* With Omron V.V.T (Virtual Valve Technology) nebulizer kit:

- Highly efficient nebulization with reduced in medication wastage
- Safe and easy to clean

- * Compact
- * Lightweight
- * Low noise

Omron V.V.T. Nebulizer kit

The valves in conventional nebulizers allow exhalation through mouthpiece. User can synchronize breathing with nebulization, hence reducing medication wastage. However, silicon valves are used in these nebulizers.

Now with Omron's unique technology, Omron nebulizers offer benefits of valve function without silicon valves.

Advantages of no silicon valves:

- Safe and easy-to-clean nebulizer kit
- No small parts which can be easily lost and swallowed by young children
- Meets hygiene standards to avoid infection (autoclavable and ease of disinfection)

Technical Specifications: -

* Warranty: 1 Year

