

JATI AUTOMASI ENTERPRISE

(Co. no 001773550-P)

No. 555 C Jalan Samudra Utara 1, Taman Samudra, 68100 Batu Caves, Selangor Darul Ehsan.
Tel : 03-61859326 Fax : 03-61871294

PNEUMATIC ROTATING STOOL MODEL WM 540

Features:

- Height adjustable by pneumatic system (gas spring)
- Base of aluminum 5 legged with 50mm swivel castors
- PVC covered cushion foam seat

Make:

- Constructed of tubular steel frame
- Finishing with a layer of chip and scratch resistant baked epoxy coat.

Size:

- Height range: 17" to 22" (430mm – 560mm)
- Seat of diameter 14" (360mm) x 3" (75mm) thick



ROTATING STOOL MODEL WM 510

Features :

- Height adjustable by revolving/rotating of the seat WITH STOPPER to prevent from rotate out
- Seat of PVC covered with PU foam
- Mounted on 4 x 50mm non-marking castors

Make :

- Constructed of tubular steel leg frame coated with hard epoxy powder coating

Size :

- Seat dimension: 14" (36cm) diameter x 3" (7.5cm) thick

Rotating Height Adjustable:

- 480-640mmH



EXAMINATION STEPS STOOL MODEL WM 101

- Single step
- Constructed of tubular steel leg frame coated with hard epoxy powder coating
- Top of plywood covered with non-slipped PVC matt.
- Mounted on rubber studs.
- Size: 508mmL x 295mmW x 170mmH

